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ESTABILIDADE

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5.3.2 Storage of specimens

Specimens should be correctly collected and delivered as quickly as possible to the laboratory. Every effort must be made to organise and expedite specimen transportation and processing. Although TB bacilli can survive in sputum for one week in the absence of preservatives, the probability of successfully culturing the bacilli decreases with time and this is especially critical for paucibacillary specimens. If specimens cannot be transported to the laboratory within one hour, it is recommended to store them at 4°C. This does not apply to whole blood specimens, which are not to be refrigerated. On arrival at the laboratory, specimens should again be refrigerated until they can be processed. The delay between collection and inoculation should not exceed seven days.

European Centre for Disease Prevention and Control. Handbook on tuberculosis laboratory diagnostic methods in the European Union – Updated 2022. Stockholm: ECDC; 2023.

delay in transport of the specimen to the laboratory was the only independent factor significantly associated with increase in culture contamination.

Ellappan, Kalaiarasan et al. “Evaluation of factors influencing Mycobacterium tuberculosis complex recovery and contamination rates in MGIT960.” *The Indian journal of tuberculosis* vol. 67,4 (2020): 466-471. doi:10.1016/j.ijtb.2020.07.016


UN 3373 “Biological Substance, Category B”

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TRANSPORTE

Packing and labelling of Category B infectious substances

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
1 Carbol Fuchsin 1 000x	2 Report	3 Fluorechent microscopy magnification*		
		250x	450x	630x
0	No acid-fast bacilli seen	0	0	0
1-9 / 100 fields	Report exact count	Divide observed count by 10	Divide observed count by 4	Divide observed count by 2
10-99/100 fields	+1			
1-10/field	+2			
>10 / field	+3			

* To adjust for altered magnification of fluorescent microscope divide the number of organisms seen by the factor provided and refer to column 1 for range and column 2 for what to report.

World Health Organization. Laboratory Service in Tuberculosis Control. Part II:Microscopy. Geneva,1998

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


MÉTODOS MOLECULARES - TAAN

In 2008, WHO approved the use of commercial LPAs for detecting MTBC in combination with resistance to rifampicin and isoniazid in sputum smear-positive specimens (direct testing) and in cultured isolates of MTBC (indirect testing).

WHO consolidated guidelines on tuberculosis. Module 3: diagnosis - rapid diagnostics for tuberculosis detection, 2021 update

Rapid Implementation of the Xpert MTB/RIF diagnostic test

Technical and operational 'How-to'
Practical considerations

March 2011

- The Xpert MTB/RIF assay provides the foundation for a revolution in the diagnosis of TB and drug-resistant TB.
- Global roll-out of the Xpert MTB/RIF assay and associated GeneXpert instruments has started. By 30 June 2011, 26 of the 145 countries that are eligible to purchase instruments and Xpert MTB/RIF cartridges at concessional prices had done so.

Global tuberculosis control: WHO report 2011 (WHO/HTM/TB/2011.16)

MÉTODOS MOLECULARES - TAAN


WHO Meeting Report of a Technical Expert
 Consultation: Non-inferiority analysis of Xpert MTB/RIF
 Ultra compared to Xpert MTB/RIF

2017

This has resulted in Ultra having a limit of detection (LOD) of 16 bacterial colony forming units (cfu) per ml (compared to 114 cfu per ml for Xpert MTB/RIF).

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
Xpert® MTB/XDR




Table 2.3.1. Drug-related gene targets, codon regions and nucleotide sequences that determine presence of variants associated with drug resistance in the Xpert MTB/XDR assay (35)

Drug	Gene target	Codon regions	Nucleotide
Isoniazid	<i>inhA</i> promoter	Not applicable	-1 to -32 intergenic region
	<i>katG</i>	311-319	939-957
	<i>fabG1</i>	199-210	597-630
	<i>oxyR-ahpC</i> intergenic region	Not applicable	-5 to -50 intergenic region (or -47 to -92 ^a)
Ethionamide	<i>inhA</i> promoter	Not applicable	-1 to -32 intergenic region
Fluoroquinolones	<i>gyrA</i>	87-95	261-285
	<i>gyrB</i>	531-544 (or 493-505 ^a)	1596-1632
Amikacin, capreomycin and kanamycin	<i>rrs</i>	Not applicable	1396-1417
Amikacin and kanamycin	<i>eis</i> promoter	Not applicable	-6 to -42 intergenic region

^a Codon numbering system according to Camus et al. (2002) (36).

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

fluoroquinolonas

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MÉTODOS MOLECULARES - TAAN

It is important to note that, like other molecular tests currently endorsed by WHO, LPAs have some limitations:

- Although LPAs can detect the mutations most frequently identified in resistant strains, some mutations that confer resistance are outside the regions covered by the test and resistance cannot be completely excluded, even in the presence of all WT probes. Thus, in some cases additional phenotypic DST may be necessary for a full assessment of the presence of a resistant strain.
- LPA is less efficient than conventional culture-based DST in detecting resistance in samples that harbour both drug-susceptible and -resistant bacteria (i.e., heteroresistance). Specifically, LPA can be used to identify resistant bacteria with mutations detected by the MUT probes if resistant bacteria represent at least 5% of the total population; however, resistant bacteria with mutations inferred by the absence of WT probes would probably be missed if the resistant population represents less than 95% of the total bacterial population (13, 14).

Line probe assays for detection of drug-resistant tuberculosis: interpretation and reporting manual for laboratory staff and clinicians


ISBN 978-92-4-004666-5 (electronic version)

ISBN 978-92-4-004667-2 (print version)

© World Health Organization 2022

Diagnostic tests for TB disease have improved substantially in recent years. There are now several rapid molecular tests recommended by WHO as the initial diagnostic test for TB, some of which can detect drug resistance simultaneously (3). These tests can be used at the lower levels of the health system. A point-of-care lateral-flow test performed on urine is also recommended by WHO; its main use is to assist with diagnosis of TB in people with advanced HIV disease, in combination with rapid molecular tests. There are additional rapid molecular tests specifically for the detection of resistance to a variety of first- and second-line anti-TB drugs, while sequencing technologies can be used to provide a comprehensive individual profile of drug resistance. The older method of sputum smear microscopy (developed >100 years ago) is still used for TB diagnosis in low and middle-income countries but is increasingly being replaced with rapid tests.

Culture testing remains the reference standard for TB diagnosis. In addition, culture is required for the detection of resistance to newer anti-TB drugs and may also be used as a confirmatory test in settings and situations in which people have a low pre-test probability of having TB disease. Following diagnosis, culture or smear (as opposed to rapid molecular tests) are necessary to monitor an individual's response to treatment.

Global tuberculosis report 2025. Geneva: World Health Organization; 2025. Licence: CC BY-NC-SA 3.0 IGO.

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

- método lento

EXAME CULTURAL

Figure 5.1 Minimal examination schedule

European Centre for Disease Prevention and Control. Handbook on TB laboratory diagnostic methods for the European Union, Stockholm: ECDC; 2016

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


LF - LAM

Tests based on the detection of mycobacterial lipoarabinomannan (LAM) antigen in urine have emerged as potential point-of-care tests for TB. The currently available urinary LAM assays have suboptimal sensitivity, and are therefore not suitable as diagnostic tests for TB in all populations. However, unlike traditional diagnostic methods, urinary LAM assays demonstrate improved sensitivity for the diagnosis of TB among individuals coinfecting with HIV. The estimated sensitivity is even greater in patients with lower CD4 cell counts. Currently, the urine lateral flow LAM assay (LF-LAM) strip-test – the Alere Determine™ TB LAM Ag, United States of America (USA), hereafter referred to as AlereLAM – is the only commercially available urinary LAM test.

Remarks

All patients with signs and symptoms of pulmonary TB who are capable of producing sputum should have as their initial diagnostic test at least one sputum specimen submitted for Xpert® MTB/RIF (Ultra) assay. This also includes children and adolescents living with HIV who are able to provide a sputum sample.

LF-LAM should be used as an add-on to clinical judgement in combination with other tests. It should not be used as a replacement or triage test.

LF - LAM

WHO policy recommendations

In inpatient settings, WHO strongly recommends using LF-LAM to assist in the diagnosis of active TB in HIV-positive adults, adolescents and children:

- with signs and symptoms of TB (pulmonary and/or extrapulmonary) (**strong** recommendation; moderate certainty in the evidence about the intervention effects) [1];³ or
- with advanced HIV disease⁴ or who are seriously ill⁵ (**strong** recommendation; moderate certainty in the evidence about the intervention effects) [1]; or
- irrespective of signs and symptoms of TB and with a CD4 cell count of less than 200 cells/mm³ (**strong** recommendation; moderate certainty in the evidence about the intervention effects) [2].

In outpatient settings, WHO suggests using LF-LAM to assist in the diagnosis of active TB in HIV-positive adults, adolescents and children:

- with signs and symptoms of TB (pulmonary and/or extrapulmonary) or seriously ill (**conditional** recommendation; low certainty in the evidence about test accuracy) [3]; and
- irrespective of signs and symptoms of TB and with a CD4 cell count of less than 100 cells/mm³ (**conditional** recommendation; very low certainty in the evidence about test accuracy) [4].

LF - LAM

In outpatient settings, WHO recommends **against** using LF-LAM to assist in the diagnosis of active TB in HIV-positive adults, adolescents and children:

- without assessing TB symptoms (**strong** recommendation; very low certainty in the evidence about test accuracy) [5];
- without TB symptoms and unknown CD4 cell count or without TB symptoms and CD4 cell count greater than or equal to 200 cells/mm³ (**strong** recommendation; very low certainty in the evidence about test accuracy) [6]; and
- without TB symptoms and with a CD4 cell count of 100–200 cells/mm³ (**conditional** recommendation; very low certainty in the evidence about test accuracy) [7].

Fig. 4. Procedimento do ensaio Alere Determine TB LAM Ag

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


PCR (MTC e MNT)

WHO Meeting Report of a Technical Expert Consultation: Non-inferiority analysis of Xpert MTB/RIF Ultra compared to Xpert MTB/RIF

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





Tuberculosis (TB) is a preventable and usually curable disease. Nonetheless, more than 10 million people continue to fall ill with TB every year and more than 1 million die from the disease, making it the world's leading cause of death from a single infectious agent and among the top 10 causes of death worldwide.

Yes! We can End TB!

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anabela.silva@insa.min-saude.pt